


NANOPASS®

33

The World's Thinnest & First Double-Tapered Pen Needle



NANOPASS® 33 won the „Good Design Grand Prize“ (2005) of the Japan Industrial Design Promotion Organization (JIDPO)

NANOPASS®

33



Product code:	TN*3305M
Needle tip Ø:	33G (0.2mm)
Needle length:	5mm
Packaging size:	100 pieces

"TYPE A" ISO 11608-2
 CE 0197  LATEX



Any questions?
www.terumo-diabetes.com
toll free number 0800 988 4617

www.terumo-europe.com

DI 89GB-0410FD-11 (08.10)E

NANOPASS®

33

Needle for Pen-injectors

The World's Thinnest Pen Needle



TERUMO®
We keep life flowing



The advantages

- Double Tapered**
new technology
optimal insulin flow
- 33G (Ø 0.2mm)**
33% thinner than
conventional 30G
needle
- Less Ouch**
easy injection
effective therapy
- One For All**
compatible with all
pen injectors (Type A)


World's Thinnest Pen Needle 33G (Ø 0.2mm)

A thinner pen needle causes less pain. The diameter of the tip of NANOPASS® 33 is only 0.2mm.

That means up to 33% thinner than a conventional 30G needle.


NANOPASS® 33 fits neatly into your personal life

- less pain
- easy injection
- optimal insulin flow



First Double-Tapered Pen Needle

The thinner a needle is, the greater its injection resistance. The double-tapered structure of NANOPASS® 33 has gone beyond conventional technology.



NANOPASS® 33 reduces the insulin injection force noticeably.

Try it and see for yourself!

